

Pneumonia Detection Using Chest X-Ray with Deep Learning Project Setup

1. Hardware Requirements:

- Memory -4.0 GB
- Processor Type -64 bit PROCESSOR
- Operating system -windows 8/windows 10

2. Software Requirements:

Keras



Jupyter Notebook



TensorFlow



Numpy



Pandas



OpenCV



Pickle



Tkinter



3. Dataset:

- Download the dataset from [Kaggle](https://www.kaggle.com)

4. Execution:

- Move to project directory using change directory command.
For Example: C:\Users\ELCOT\Desktop>cd pneumonia_detection
- Now run the command **jupyter notebook <notebook name>** in command prompt.

For Example: C:\Users\ELCOT\Desktop\pneumonia_detection>jupyter notebook Pneumonia_Detection_project.ipynb

5. Program Code:

LIBRARY FILES

```
import tensorflow as tf
import numpy as np
import os
import sys
import cv2
```

```
import matplotlib.pyplot as plt
import pickle
import random
import pandas as pd
```

```
import tensorflow as tf
from tensorflow.keras import Sequential
from tensorflow.keras.layers import Dense,Dropout,Activation,
Flatten,Conv2D,MaxPooling2D
import pickle
```

```
from sklearn.model_selection import train_test_split
```

```
import matplotlib.pyplot as plt
%matplotlib inline
```

PREPROCESSING

```
class MasterImage(object):
```

```
    def __init__(self,PATH="", IMAGE_SIZE = 50, CATEGORIES=[]):
        self.PATH = PATH
        self.IMAGE_SIZE = IMAGE_SIZE
        self.CATEGORIES = CATEGORIES
        self.image_data = []
        self.x_data = []
        self.y_data = []
```

```
    def Process_Image(self):
```

```
        """
        Return Numpy array of image
        :return: X_Data, Y_Data
        """
```

```
        for categories in self.CATEGORIES:                                # Iterate
over categories
```

```
            train_folder_path = os.path.join(self.PATH, categories)      #
Folder Path
```

```
            class_index = self.CATEGORIES.index(categories)              # this
will get index for classification
```

```
            for img in os.listdir(train_folder_path):                    # This will
iterate in the Folder
```

```
                new_path = os.path.join(train_folder_path, img)          # image
Path
```

```
                try:              # if any image is corrupted
```

```

        image_data_temp =
cv2.imread(new_path,cv2.IMREAD_GRAYSCALE)           # Read Image as
numbers
        image_temp_resize =
cv2.resize(image_data_temp,(self.IMAGE_SIZE,self.IMAGE_SIZE))
        self.image_data.append([image_temp_resize,class_index])
    except:
        pass

data = np.asanyarray(self.image_data)

# Iterate over the Data
for x in data:
    self.x_data.append(x[0])    # Get the X_Data
    self.y_data.append(x[1])    # get the label

X_Data = np.asarray(self.x_data) / (255.0)    # Normalize Data
Y_Data = np.asarray(self.y_data)

return X_Data,Y_Data

def pickle_image(self):

    """
    :return: None Creates a Pickle Object of DataSet
    """

    X_Data,Y_Data = self.Process_Image()

    pickle_out = open('X_Data','wb')
    pickle.dump(X_Data, pickle_out)
    pickle_out.close()

    pickle_out = open('Y_Data', 'wb')
    pickle.dump(Y_Data, pickle_out)
    pickle_out.close()

    print("Pickled Image Successfully ")

    return X_Data,Y_Data

def load_dataset(self):

    try:
        X_Temp = open('X_Data','rb')
        X_Data = pickle.load(X_Temp)

        Y_Temp = open('Y_Data','rb')
        Y_Data = pickle.load(Y_Temp)

```

```

        print('Reading Dataset from Pickle Object')

        return X_Data,Y_Data

    except:
        print('Could not Found Pickle File ')
        print('Loading File and Dataset .....')
        X_Data,Y_Data = self.pickle_image()
        return X_Data,Y_Data
a=MasterImage(PATH='/Users/ELCOT/Desktop/prjct/chest_xray/chest_xray/train',
               IMAGE_SIZE=80,
               CATEGORIES=['NORMAL','PNEUMONIA'])
X_Data,Y_Data = a.load_dataset()
X_Data = X_Data.reshape(-1,80,80,1)
print(X_Data.shape)
SPLITTING THE TRAINING AND TESTING DATA
X_Train, X_Test, Y_Train, Y_Test = train_test_split(X_Data, Y_Data,
test_size=0.3,random_state=101)
TRAINING THE MODEL
model = Sequential()
model.add(Conv2D(200, (3, 3), input_shape=X_Data.shape[1:]))
model.add(Activation('relu'))
model.add(MaxPooling2D(pool_size=(2, 2)))
model.add(Conv2D(100, (3, 3)))
model.add(Activation('relu'))
model.add(MaxPooling2D(pool_size=(2, 2)))
model.add(Flatten())
model.add(Dense(80))
model.add(Dense(1))
model.add(Activation('sigmoid'))
model.compile(loss='binary_crossentropy',
              optimizer='adam',
              metrics=['accuracy'])
model.fit(X_Data, Y_Data, batch_size=40, epochs=5, validation_split=0.3)
model.save('model.h5')
EVALUATING THE MODEL
model.evaluate(X_Test,Y_Test,batch_size=40)
LOADING THE MODEL
model = tf.keras.models.load_model('model.h5')
USER INTERFACE
import tkinter as tk
from tkinter.filedialog import askopenfilename
import shutil
import os
import sys
from PIL import Image, ImageTk

window = tk.Tk()

```

```

window.title("PNEUMONIA DETECTION")

window.geometry("500x510")
window.configure(background="black")

title = tk.Label(text="Click below to choose picture for testing disease....",
background="black", fg="white", font=("", 15))
title.grid()
def exit():
    window1.destroy()

def prepare(filepath):
    training_date = []

    img_array = cv2.imread(filepath,cv2.IMREAD_GRAYSCALE)
    new_array = cv2.resize(img_array,(80,80))
    new_image = new_array.reshape(-1,80,80,1)
    return new_image

def analysis():

    #model = tf.keras.models.load_model('model.h5')
    test = model.predict([prepare(filepath=filepath)])
    CATEGORIES=['NORMAL','PNEUMONIA']
    print(CATEGORIES[int(test[0][0])])
    #model_out = model.predict([data])[0]
    message = tk.Label(text='Status: '+CATEGORIES[int(test[0][0])],
background="black",
fg="white", font=("", 15))
    message.grid(column=0, row=3, padx=10, pady=10)
    button = tk.Button(text="Exit", command=exit)
    button.grid(column=0, row=9, padx=20, pady=20)

def openphoto():
    global filepath
    global window1
    window.destroy()
    window1 = tk.Tk()

    window1.title("PNEUMONIA DTECTION")

    window1.geometry("500x520")
    window1.configure(background="black")

    # C:/Users/sagpa/Downloads/images is the location of the image which you want
    to test..... you can change it according to the image location you have

```


```

filepath =
askopenfilename(initialdir='C:/Users/ELCOT/Desktop/prjct/chest_xray', title='Select
image for analysis ',
                filetypes=[('image files', '.jpeg')])
load = Image.open(filepath)
load= load.resize((250,490),Image.ANTIALIAS)
render = ImageTk.PhotoImage(load)
img = tk.Label(image=render, height="250", width="490")
img.image = render
img.place(x=0, y=0)
img.grid(column=0, row=0, padx=5, pady = 5)
button2 = tk.Button(text="Analyse Image", command=analysis)
button2.grid(column=0, row=2, padx=10, pady = 10)
window1.mainloop()
button1 = tk.Button(text="Get Photo", command = openphoto)
button1.grid(column=0, row=1, padx=10, pady = 10)
window.mainloop()


```

6. Result:

- **Home Page**



- **Output For Normal CXR:**



- **Output For Pneumonia CXR**

